

For Research Use Only

PD-L2/B7-DC Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68795-1-PBS

Basic Information

Catalog Number:

68795-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC074766

GeneID (NCBI):

80380

UNIPROT ID:

Q9BQ51

Full Name:

programmed cell death 1 ligand 2

Calculated MW:

273 aa, 31 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H9C10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50162-1, PD-L2/B7-DC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68795-1-PBS. Detection antibody: 68795-2-PBS. Standard:Eg0507. Range: 0.781-100 ng/mL