

For Research Use Only

# MPO Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68797-2-PBS

## Basic Information

Catalog Number:

68797-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC130476

GeneID (NCBI):

4353

UNIPROT ID:

P05164

Full Name:

myeloperoxidase

Calculated MW:

745 aa, 84 kDa

Purification Method:

Protein G purification

CloneNo.:

1D7B10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50164-1, MPO Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68797-1-PBS. Detection antibody: 68797-2-PBS. Standard:Eg0511. Range: 0.781-100 ng/mL.