

For Research Use Only

SCN9A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68798-1-PBS

Basic Information

Catalog Number:

68798-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33724

GenBank Accession Number:

XM_011511616.2

GeneID (NCBI):

6335

UNIPROT ID:

Q15858

Full Name:

sodium channel, voltage-gated, type

IX, alpha subunit

Calculated MW:

226

Purification Method:

Protein G Magarose purification

CloneNo.:

3E4D4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50165-1, SCN9A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68798-1-PBS. Detection antibody: 68798-2-PBS. Standard:Ag33724. Range: 0.781-100 ng/mL.