

For Research Use Only

# PARK7/DJ-1 Monoclonal antibody

Catalog Number: 68802-1-Ig



## Basic Information

Catalog Number:

68802-1-Ig

Size:

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2287

GenBank Accession Number:

BC008188

GeneID (NCBI):

11315

UNIPROT ID:

Q99497

Full Name:

Parkinson disease (autosomal recessive, early onset) 7

Calculated MW:

189 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

1E5C3

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

