

For Research Use Only

KCNMB3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68804-1-PBS

Basic Information

Catalog Number:

68804-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31765

GenBank Accession Number:

NM_171828

GeneID (NCBI):

27094

UNIPROT ID:

Q9NPA1

Full Name:

potassium large conductance

calcium-activated channel, subfamily

M beta member 3

Calculated MW:

32KD

Purification Method:

Protein G Magarose purification

CloneNo.:

3B2D1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50171-1, KCNMB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68804-1-PBS. Detection antibody: 68804-2-PBS. Standard:Ag31765. Range: 0.16-20 ng/mL.