

For Research Use Only

IL-12/IL-23 p40 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68805-2-PBS



Basic Information

Catalog Number:

68805-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC067498

GeneID (NCBI):

3593

UNIPROT ID:

P29460

Full Name:

interleukin 12B (natural killer cell
stimulatory factor 2, cytotoxic
lymphocyte maturation factor 2, p40)

Calculated MW:

328 aa, 37 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F9E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50172-1, IL12B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68805-1-PBS. Detection antibody: 68805-2-PBS. Standard:Eg0431. Range: 1.563-100 ng/mL.