

For Research Use Only

# HHIP Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68807-2-PBS

## Basic Information

Catalog Number:

68807-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2241

GenBank Accession Number:

BC025311

GeneID (NCBI):

64399

UNIPROT ID:

Q96QV1

Full Name:

hedgehog interacting protein

Calculated MW:

700 aa, 79 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3F2F5

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50174-1, HHIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68807-1-PBS. Detection antibody: 68807-2-PBS. Standard:Ag2241. Range: 1.563-100 ng/mL.