

For Research Use Only

Vinculin Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68809-2-PBS

Basic Information

Catalog Number:

68809-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG3705

GenBank Accession Number:

BC039174

GeneID (NCBI):

7414

UNIPROT ID:

P18206

Full Name:

vinculin

Calculated MW:

1133 aa, 124 kDa

Purification Method:

Protein A purification

CloneNo.:

2A11D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50176-1, Vinculin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68809-1-PBS. Detection antibody: 68809-2-PBS. Standard:Ag3705. Range: 0.098-6.25 ng/mL