

For Research Use Only

ATF1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68813-1-PBS

Basic Information

Catalog Number: 68813-1-PBS	GenBank Accession Number: BC029619	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 466	CloneNo.: 2H5C9
Source: Mouse	UNIPROT ID: P18846	
Isotype: IgG1	Full Name: activating transcription factor 1	
Immunogen Catalog Number: AG34935	Calculated MW: 271 aa, 29 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50187-1, ATF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68813-1-PBS. Detection antibody: 68813-2-PBS. Standard:Ag34935. Range: 0.098-12.5 ng/mL