

For Research Use Only

ATF1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68813-2-PBS

Basic Information

Catalog Number:

68813-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG34935

GenBank Accession Number:

BC029619

GeneID (NCBI):

466

UNIPROT ID:

P18846

Full Name:

activating transcription factor 1

Calculated MW:

271 aa, 29 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D9C7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50187-1, ATF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68813-1-PBS. Detection antibody: 68813-2-PBS. Standard:Ag34935. Range: 0.098-12.5 ng/mL