

For Research Use Only

CPA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68826-2-PBS



Basic Information

Catalog Number:

68826-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8579

GenBank Accession Number:

BC005279

GeneID (NCBI):

1357

UNIPROT ID:

P15085

Full Name:

carboxypeptidase A1 (pancreatic)

Calculated MW:

419 aa, 47 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E9D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Cytometric bead array sample test of MP51787-1, CPA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68826-2-PBS. Detection antibody: 68826-3-PBS.



Cytometric bead array sample test of MP51787-2, CPA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68826-2-PBS. Detection antibody: 68826-4-PBS.



Cytometric bead array standard curve of MP51787-1, CPA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68826-2-PBS. Detection antibody: 68826-3-PBS. Standard:Ag8579. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP51787-2, CPA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68826-2-PBS. Detection antibody: 68826-4-PBS. Standard:Ag8579. Range: 0.781-100 ng/mL