

For Research Use Only

# ADAMTS12 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68829-4-PBS

## Basic Information

Catalog Number:

68829-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34544

GenBank Accession Number:

BC058841

GeneID (NCBI):

81792

UNIPROT ID:

P58397

Full Name:

ADAM metalloproteinase with  
thrombospondin type 1 motif, 12

Calculated MW:

178 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B9E1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50206-2, ADAMTS12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68829-3-PBS. Detection antibody: 68829-4-PBS. Standard:Ag34544. Range: 0.098-25 ng/mL.