

For Research Use Only

TREM2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68831-1-PBS

Basic Information

Catalog Number:

68831-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC032362

GeneID (NCBI):

54209

UNIPROT ID:

Q9NZC2

Full Name:

triggering receptor expressed on myeloid cells 2

Calculated MW:

222 aa, 25 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C9G3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50208-1, TREM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68831-1-PBS. Detection antibody: HRP-conjugated 68831-2-PBS. Standard:Eg0747. Range: 0.049-50 ng/mL.