

For Research Use Only

# CALM1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68833-2-PBS

## Basic Information

Catalog Number:

68833-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG35001

GenBank Accession Number:

BC000454

GeneID (NCBI):

801

UNIPROT ID:

PODP23

Full Name:

calmodulin 1 (phosphorylase kinase, delta)

Calculated MW:

17 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G9D7

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50211-1, CALM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68833-1-PBS. Detection antibody: HRP-conjugated 68833-2-PBS. Standard:Ag35001. Range: 0.049-50 ng/mL.