

For Research Use Only

# IGFBP4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68834-2-PBS

## Basic Information

Catalog Number:

68834-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG34800

GenBank Accession Number:

BC016041

GeneID (NCBI):

3487

UNIPROT ID:

P22692

Full Name:

insulin-like growth factor binding protein 4

Calculated MW:

28 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3B10E2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50212-1, IGF1BP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68834-1-PBS. Detection antibody: HRP-conjugated 68834-2-PBS. Standard:Ag34800. Range: 0.049-50 ng/mL.