

For Research Use Only

IGFBP4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68834-4-PBS

Basic Information

Catalog Number:

68834-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG34800

GenBank Accession Number:

BC016041

GeneID (NCBI):

3487

UNIPROT ID:

P22692

Full Name:

insulin-like growth factor binding protein 4

Calculated MW:

28 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A11F5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50212-2, IGFBP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68834-3-PBS. Detection antibody: HRP-conjugated 68834-4-PBS. Standard: Ag34800. Range: 0.049-50 ng/mL.