

For Research Use Only

ICOS Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68836-2-PBS

Basic Information

Catalog Number: 68836-2-PBS	GenBank Accession Number: BC028210	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 29851	CloneNo.: 3D10E12
Source: Mouse	UNIPROT ID: Q9Y6W8	
Isotype: IgG1	Full Name: inducible T-cell co-stimulator	
Immunogen Catalog Number: AG27722	Calculated MW: 199 aa, 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50213-1, ICOS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68836-1-PBS. Detection antibody: 68836-2-PBS. Standard:Ag27722. Range: 0.098-25 ng/mL