

For Research Use Only

FOXP3 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68838-1-PBS

Basic Information

Catalog Number:

68838-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17571

GenBank Accession Number:

BC113401

GeneID (NCBI):

50943

UNIPROT ID:

Q9BZS1

Full Name:

forkhead box P3

Calculated MW:

431 aa, 47 kDa

Purification Method:

Protein G purification

CloneNo.:

1H6F3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50216-1, FOXP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68803-2-PBS. Detection antibody: 68838-1-PBS. Standard:Ag17571. Range: 0.781-100 ng/mL.