

For Research Use Only

Neurofibromin 2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68839-2-PBS



Basic Information

Catalog Number: 68839-2-PBS	GenBank Accession Number: BC020257	Purification Method: Protein G purification
Size: 1000 ug/ml	GeneID (NCBI): 4771	CloneNo.: 2B10H6
Source: Mouse	UNIPROT ID: P35240	
Isotype: IgG1	Full Name: neurofibromin 2 (merlin)	
Immunogen Catalog Number: AG16265	Calculated MW: 595 aa, 70 kDa	

Applications

Tested Applications:
IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50217-1, NF2/Neurofibromin 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68839-1-PBS. Detection antibody: 68839-2-PBS. Standard:Ag16265. Range: 0.098-100 ng/mL