

For Research Use Only

# Neurofibromin 2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68839-6-PBS



## Basic Information

Catalog Number:

68839-6-PBS

Size:

1000 ug/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16265

GenBank Accession Number:

BC020257

GeneID (NCBI):

4771

UNIPROT ID:

P35240

Full Name:

neurofibromin 2 (merlin)

Calculated MW:

595 aa, 70 kDa

Purification Method:

Protein G purification

CloneNo.:

2G10E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50217-3, NF2/Neurofibromin 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68839-5-PBS. Detection antibody: 68839-6-PBS. Standard:Ag16265. Range: 0.391-100 ng/mL