

For Research Use Only

BAP1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68840-2-PBS

Basic Information

Catalog Number:

68840-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27333

GenBank Accession Number:

BC001596

GeneID (NCBI):

8314

UNIPROT ID:

Q92560

Full Name:

BRCA1 associated protein-1 (ubiquitin
carboxy-terminal hydrolase)

Calculated MW:

81-91 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C8E3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50223-1, BAP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68840-1-PBS. Detection antibody: 68840-2-PBS. Standard:Ag27333. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50223-2, BAP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68840-3-PBS. Detection antibody: 68840-2-PBS. Standard:Ag27333. Range: 0.098-25 ng/mL.