

For Research Use Only

# PROS1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68843-1-PBS



## Basic Information

Catalog Number: 68843-1-PBS	GenBank Accession Number: BC015801	Purification Method: Protein G purification
Concentration: 1000 µg/ml	GeneID (NCBI): 5627	CloneNo.: 1B2E3
Source: Mouse	UNIPROT ID: P07225	
Isotype: IgG1	Full Name: protein S (alpha)	
Immunogen Catalog Number: AG10503	Calculated MW: 75 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50227-1, PROS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68843-1-PBS. Detection antibody: 68843-2-PBS. Standard:Ag10503. Range: 1.563-100 ng/mL.