

For Research Use Only

ZC3H12D Monoclonal antibody, PBS Only



Catalog Number: 68847-1-PBS

Basic Information

Catalog Number: 68847-1-PBS	GenBank Accession Number: BC157832	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 340152	CloneNo.: 1E2A2
Source: Mouse	UNIPROT ID: A2A288	
Isotype: IgG2b	Full Name: zinc finger CCCH-type containing 12D	
Immunogen Catalog Number: AG21770	Calculated MW: 527 aa, 58 kDa	

Applications

Tested Applications:
ELISA, WB

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68847-1-Ig (ZC3H12D antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68847-1-PBS in a different storage buffer formulation.