

For Research Use Only

# PINK1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68848-2-PBS



## Basic Information

Catalog Number: 68848-2-PBS	GenBank Accession Number: BC028215	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 65018	CloneNo.: 2F9G12
Source: Mouse	UNIPROT ID: Q9BXM7	
Isotype: IgG2b	Full Name: PTEN induced putative kinase 1	
Immunogen Catalog Number: AG19825	Calculated MW: 581 aa, 63 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50236-1, PINK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68848-1-PBS. Detection antibody: 68848-2-PBS. Standard:Ag19825. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50236-2, PINK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68848-3-PBS. Detection antibody: 68848-2-PBS. Standard:Ag19825. Range: 0.391-100 ng/mL.