

For Research Use Only

Chromogranin A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68851-1-PBS



Basic Information

Catalog Number:

68851-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19780

GenBank Accession Number:

BC006459

GeneID (NCBI):

1113

UNIPROT ID:

P10645

Full Name:

chromogranin A (parathyroid
secretory protein 1)

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

1D11G12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50241-1, CHGA /Chromogranin A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68851-1-PBS. Detection antibody: 60135-3-PBS. Standard:Ag19780. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50241-3, CHGA/Chromogranin A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68851-1-PBS. Detection antibody: 60135-2-PBS. Standard:Ag19780. Range: 1.563-100 ng/mL.