

For Research Use Only

# RPL22 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68854-1-PBS

## Basic Information

|                                      |                                       |  |
|--------------------------------------|---------------------------------------|--|
| Catalog Number:<br>68854-1-PBS       | GenBank Accession Number:<br>BC058887 | Purification Method:<br>Protein G purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>6146                | CloneNo.:<br>2A3C10                            |
| Source:<br>Mouse                     | UNIPROT ID:<br>P35268                 |  |
| Isotype:<br>IgG1                     | Full Name:<br>ribosomal protein L22   |  |
| Immunogen Catalog Number:<br>AG30115 | Calculated MW:<br>15 kDa              |  |

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50247-1, RPL22 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68854-1-PBS. Detection antibody: 68854-2-PBS. Standard:Ag30115. Range: 1.563-50 ng/mL