

For Research Use Only

GS15 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68856-1-PBS

Basic Information

Catalog Number:

68856-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5506

GenBank Accession Number:

BC032779

GeneID (NCBI):

51272

UNIPROT ID:

Q9NYM9

Full Name:

blocked early in transport 1 homolog
(S. cerevisiae)-like

Calculated MW:

102 aa, 11 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G3B11

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50251-1, GS15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68856-1-PBS. Detection antibody: 68856-2-PBS. Standard:Ag5506. Range: 0.195-50 ng/mL



Cytometric bead array standard curve of MP50251-2, GS15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68856-1-PBS. Detection antibody: 68856-3-PBS. Standard:Ag5506. Range: 0.098-6.25 ng/mL