

For Research Use Only

# IRX4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68860-2-PBS

## Basic Information

Catalog Number:

68860-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12205

GenBank Accession Number:

BC110912

GeneID (NCBI):

50805

UNIPROT ID:

P78413

Full Name:

iroquois homeobox 4

Calculated MW:

519 aa, 54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B5F4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50254-1, IRX4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68860-1-PBS. Detection antibody: 68860-2-PBS. Standard:Ag12205. Range: 0.391-25 ng/mL



Cytometric bead array standard curve of MP50254-2, IRX4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68860-3-PBS. Detection antibody: 68860-2-PBS. Standard:null. Range: 0.391-100 ng/mL.