

For Research Use Only

FGFBasic-TS Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68867-1-PBS

Basic Information

Catalog Number:

68867-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1285

GenBank Accession Number:

GeneID (NCBI):

2247

UNIPROT ID:

P09038

Full Name:

fibroblast growth factor 2 (basic)

Purification Method:

Protein G purification

CloneNo.:

1B2B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50262-1, FGF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68867-1-PBS. Detection antibody: 68870-1-PBS. Standard:HZ-1285. Range: 0.098-100 ng/mL.