

For Research Use Only

MYO18A Monoclonal antibody, PBS Only



Catalog Number: 68868-1-PBS

Basic Information

Catalog Number:

68868-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6185

GenBank Accession Number:

BC039612

GeneID (NCBI):

399687

UNIPROT ID:

Q92614

Full Name:

myosin XVIII A

Calculated MW:

233 kDa

Purification Method:

Protein A purification

CloneNo.:

1A2B12

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68868-1-Ig (MYO18A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68868-1-PBS in a different storage buffer formulation.