

For Research Use Only

APRIL, TNFSF13 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68876-2-PBS



Basic Information

Catalog Number:

68876-2-PBS

Concentration:

1000 μ g/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG26281

GenBank Accession Number:

BC008042

GeneID (NCBI):

8741

UNIPROT ID:

O75888

Full Name:

tumor necrosis factor (ligand)
superfamily, member 13

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

2E10B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array of Human saliva using MP50273-1, APRIL,TNFSF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68876-1-PBS. Detection antibody: 68876-2-PBS. Saliva of two individual healthy human donors was analyzed using MP50273-1.



Cytometric bead array standard curve of MP50273-1, APRIL,TNFSF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68876-1-PBS. Detection antibody: 68876-2-PBS. Standard:Ag26281. Range: 1.563-100 ng/mL.