

For Research Use Only

# DLL4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68882-1-PBS



## Basic Information

Catalog Number:

68882-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17483

GenBank Accession Number:

BC106950

GeneID (NCBI):

54567

UNIPROT ID:

Q9NR61

Full Name:

delta-like 4 (Drosophila)

Calculated MW:

685 aa, 75 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A3H2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50281-1, DLL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68882-1-PBS. Detection antibody: 68882-2-PBS. Standard:Ag16140. Range: 6.25-100 ng/mL



Cytometric bead array standard curve of MP50281-2, DLL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68882-1-PBS. Detection antibody: 68882-3-PBS. Standard:Ag16140. Range: 0.391-12.5 ng/mL