

For Research Use Only

DLL4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68882-2-PBS



Basic Information

Catalog Number:

68882-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17483

GenBank Accession Number:

BC106950

GeneID (NCBI):

54567

UNIPROT ID:

Q9NR61

Full Name:

delta-like 4 (Drosophila)

Calculated MW:

685 aa, 75 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F11D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50281-1, DLL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68882-1-PBS. Detection antibody: 68882-2-PBS. Standard:Ag16140. Range: 6.25-100 ng/mL