

For Research Use Only

FOXA1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68896-3-PBS



Basic Information

Catalog Number:

68896-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG14502

GenBank Accession Number:

BC033890

GeneID (NCBI):

3169

UNIPROT ID:

P55317

Full Name:

forkhead box A1

Calculated MW:

473 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

1D9F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50297-2, FOXA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68896-1-PBS. Detection antibody: 68896-3-PBS. Standard:Ag14502. Range: 1.563-100 ng/mL.