

For Research Use Only

ERCC3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68901-1-PBS

Basic Information

Catalog Number:

68901-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20286

GenBank Accession Number:

BC008820

GeneID (NCBI):

2071

UNIPROT ID:

P19447

Full Name:

excision repair cross-complementing
rodent repair deficiency,
complementation group 3 (xeroderma
pigmentosum group B
complementing)

Calculated MW:

89 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G7B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50301-1, ERCC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68901-1-PBS. Detection antibody: 68901-2-PBS. Standard:Ag20286. Range: 0.781-100 ng/mL.