

For Research Use Only

ZCCHC3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68903-1-PBS

Basic Information

| | | |
|--------------------------------------|--------------------------------------------------------|------------------------------------------------|
| Catalog Number: 68903-1-PBS | GenBank Accession Number: BC069238 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 85364 | CloneNo.: 1G9G3 |
| Source: Mouse | UNIPROT ID: Q9NUD5 | |
| Isotype: IgG2a | Full Name: zinc finger, CCHC domain containing 3 | |
| Immunogen Catalog Number: AG34163 | | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50305-1, ZCCHC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68903-1-PBS. Detection antibody: 68903-2-PBS. Standard:Ag34163. Range: 0.781-200 ng/mL



Cytometric bead array standard curve of MP50305-2, ZCCHC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68903-1-PBS. Detection antibody: 68903-3-PBS. Standard:Ag34163. Range: 0.781-200 ng/mL