

For Research Use Only

ZCCHC3 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68903-2-PBS

Basic Information

Catalog Number: 68903-2-PBS	GenBank Accession Number: BC069238	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 85364	CloneNo.: 2G9B5
Source: Mouse	UNIPROT ID: Q9NUD5	
Isotype: IgG1	Full Name: zinc finger, CCHC domain containing 3	
Immunogen Catalog Number: AG34163		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50305-1, ZCCHC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68903-1-PBS. Detection antibody: 68903-2-PBS. Standard:Ag34163. Range: 0.781-200 ng/mL.