

For Research Use Only

# ZCCHC3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68903-4-PBS

## Basic Information

Catalog Number: 68903-4-PBS	GenBank Accession Number: BC069238	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 85364	CloneNo.: 1A1A2
Source: Mouse	UNIPROT ID: Q9NUD5	
Isotype: IgG2a	Full Name: zinc finger, CCHC domain containing 3	
Immunogen Catalog Number: AG34163		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50305-3, ZCCHC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68903-4-PBS. Detection antibody: 68903-5-PBS. Standard:Ag34163. Range: 6.25-200 ng/mL