

For Research Use Only

ROPN1L Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68909-3-PBS

Basic Information

Catalog Number:

68909-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34015

GenBank Accession Number:

BC014607

GeneID (NCBI):

83853

UNIPROT ID:

Q96C74

Full Name:

ropporin 1-like

Calculated MW:

230 aa, 26 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H10G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50320-2, ROPN1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68909-3-PBS. Detection antibody: 68909-2-PBS. Standard:Ag34015. Range: 0.781-100 ng/mL.