

For Research Use Only

TSLP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68911-1-PBS



Basic Information

Catalog Number:

68911-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4771

GenBank Accession Number:

BC016720

GeneID (NCBI):

85480

UNIPROT ID:

Q969D9

Full Name:

thymic stromal Lymphopoietin

Calculated MW:

15 kDa, 42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H1G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50325-1, TSLP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68911-1-PBS. Detection antibody: 68911-2-PBS. Standard:Ag4771. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50325-2, TSLP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68911-1-PBS. Detection antibody: 68911-3-PBS. Standard:Ag4771. Range: 0.098-50 ng/mL.