

For Research Use Only

WT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68912-1-PBS



Basic Information

Catalog Number:

68912-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19541

GenBank Accession Number:

BC032861

GeneID (NCBI):

7490

UNIPROT ID:

P19544

Full Name:

Wilms tumor 1

Calculated MW:

449 aa, 49 kDa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B10D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50330-1, WT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68912-1-PBS. Detection antibody: 68912-2-PBS. Standard:Ag19541. Range: 6.25-200 ng/mL



Cytometric bead array standard curve of MP50330-2, WT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68912-1-PBS. Detection antibody: 68912-3-PBS. Standard:Ag19541. Range: 6.25-200 ng/mL