

For Research Use Only

BICC1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68914-2-PBS

Basic Information

Catalog Number:

68914-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG35028

GenBank Accession Number:

BC156900

GeneID (NCBI):

80114

UNIPROT ID:

Q9H694

Full Name:

bicaudal C homolog 1 (Drosophila)

Calculated MW:

105 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2A9C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50332-1, BICC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68914-1-PBS. Detection antibody: 68914-2-PBS. Standard:Ag35028. Range: 0.391-100 ng/mL.