

For Research Use Only

# BICC1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68914-3-PBS

## Basic Information

Catalog Number:

68914-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG35028

GenBank Accession Number:

BC156900

GeneID (NCBI):

80114

UNIPROT ID:

Q9H694

Full Name:

bicaudal C homolog 1 (Drosophila)

Calculated MW:

105 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A12D10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50332-2, BICC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68914-3-PBS. Detection antibody: 68914-4-PBS. Standard:Ag35028. Range: 0.391-25 ng/mL