

For Research Use Only

BICC1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68914-4-PBS

Basic Information

Catalog Number: 68914-4-PBS	GenBank Accession Number: BC156900	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 80114	CloneNo.: 3D7B11
Source: Mouse	UNIPROT ID: Q9H694	
Isotype: IgG3	Full Name: bicaudal C homolog 1 (Drosophila)	
Immunogen Catalog Number: AG35028	Calculated MW: 105 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50332-2, BICC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68914-3-PBS. Detection antibody: 68914-4-PBS. Standard:Ag35028. Range: 0.391-25 ng/mL