

For Research Use Only

PARK7/DJ-1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68915-4-PBS

Basic Information

Catalog Number:

68915-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28526

GenBank Accession Number:

BC008188

GeneID (NCBI):

11315

UNIPROT ID:

Q99497

Full Name:

Parkinson disease (autosomal
recessive, early onset) 7

Calculated MW:

189 aa, 20 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A10G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50337-3, PARK7/DJ-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68915-4-PBS. Detection antibody: 68915-5-PBS. Standard: Ag28526. Range: 1.563-100 ng/mL.