

For Research Use Only

PARK7/DJ-1 Monoclonal antibody, PBS Only

Catalog Number: 68915-5-PBS



Basic Information

Catalog Number:

68915-5-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28526

GenBank Accession Number:

BC008188

GeneID (NCBI):

11315

UNIPROT ID:

Q99497

Full Name:

Parkinson disease (autosomal recessive, early onset) 7

Calculated MW:

189 aa, 20 kDa

Purification Method:

Protein G purification

CloneNo.:

4G9H9

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

