

For Research Use Only

Neprilysin/CD10 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68921-1-PBS



Basic Information

Catalog Number: 68921-1-PBS	GenBank Accession Number: BC106070	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4311	CloneNo.: 8D8E6
Source: Mouse	UNIPROT ID: P08473	
Isotype: IgG2b	Full Name: membrane metallo-endoropeptidase	
Immunogen Catalog Number: AG12506	Calculated MW: 80aa,9 kDa; 750aa,85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50348-1, MME,CD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60034-5-PBS. Detection antibody: 68921-1-PBS. Standard:Ag12506. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP50348-2, MME,CD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68921-2-PBS. Detection antibody: 68921-1-PBS. Standard:Ag12506. Range: 6.25-100 ng/mL