

For Research Use Only

CCDC109A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68926-3-PBS



Basic Information

Catalog Number:

68926-3-PBS

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22748

GenBank Accession Number:

BC034235

GeneID (NCBI):

90550

UNIPROT ID:

Q8NE86

Full Name:

coiled-coil domain containing 109A

Calculated MW:

35 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G11E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50349-2, CCDC109A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68926-3-PBS. Detection antibody: 68926-4-PBS. Standard:Ag22748. Range: 3.125-100 ng/mL.