

For Research Use Only

Epiregulin Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68929-1-PBS



Basic Information

Catalog Number: 68929-1-PBS	GenBank Accession Number: BC136404	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 2069	CloneNo.: 2F6G4
Source: Mouse	UNIPROT ID: O14944	
Isotype: IgG1	Full Name: epiregulin	
Immunogen Catalog Number: AG18915	Calculated MW: 169 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50352-1, EREG/Epiregulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68929-1-PBS. Detection antibody: 68929-2-PBS. Standard:Ag18915. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP50352-2, EREG/Epiregulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68929-1-PBS. Detection antibody: 68929-3-PBS. Standard:Ag18915. Range: 1.563-100 ng/mL.