

For Research Use Only

Epiregulin Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68929-2-PBS



Basic Information

Catalog Number:

68929-2-PBS

Size:

1000 μ g/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18915

GenBank Accession Number:

BC136404

GeneID (NCBI):

2069

UNIPROT ID:

O14944

Full Name:

epiregulin

Calculated MW:

169 aa, 19 kDa

Purification Method:

Protein G purification

CloneNo.:

2B4G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50352-1, EREG/Epiregulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68929-1-PBS. Detection antibody: 68929-2-PBS. Standard:Ag18915. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP50352-3, EREG /Epiregulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68929-2-PBS. Detection antibody: 68929-4-PBS. Standard:Ag18915. Range: 3.125-200 ng/mL.